

# Beta-1 adrenergic receptor

Catalog No: #AP80299



Package Size: #AP80299-1 50ug #AP80299-2 100ug #AP80299-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Beta-1 adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Beta-1 adrenoreceptor, Beta-1 adrenoceptor
Accession No.	Uniprot ID:P79148
Uniprot	P79148
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: ADRB1

Note: This product is for in vitro research use only