

Aldose reductase

Catalog No: #AP80302



Package Size: #AP80302-1 50ug #AP80302-2 100ug #AP80302-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aldose reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	AR,20-alpha-hydroxysteroid dehydrogenase,20-alpha-HSD,Aldehyde reductase
Accession No.	Uniprot ID:P16116
Uniprot	P16116
GeneID	317748;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AKR1B1

Note: This product is for in vitro research use only