

# Adenosine deaminase

Catalog No: #AP80306



Package Size: #AP80306-1 50ug #AP80306-2 100ug #AP80306-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Adenosine deaminase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Adenosine aminohydrolase,RCJMB04_9m8
Accession No.	Uniprot ID:Q5ZKP6
Uniprot	Q5ZKP6
GeneID	419194;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ADA

Note: This product is for in vitro research use only