

Protein AMBP

Catalog No: #AP80307



Package Size: #AP80307-1 50ug #AP80307-2 100ug #AP80307-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein AMBP
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Itil
Accession No.	Uniprot ID:Q07456
Uniprot	Q07456
GeneID	11699;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ambp

Note: This product is for in vitro research use only