

# ADP/ATP translocase 2

Catalog No: #AP80309



Package Size: #AP80309-1 50ug #AP80309-2 100ug #AP80309-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ADP/ATP translocase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ADP,ATP carrier protein 2,Adenine nucleotide translocator 2,ANT 2,Solute carrier family 25 member 5,Ant2
Accession No.	Uniprot ID:Q09073
Uniprot	Q09073
GeneID	25176;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Slc25a5

Note: This product is for in vitro research use only