

Disintegrin and metalloproteinase domain-containing protein 10



Catalog No: #AP80311

Orders: order@signalwayantibody.com

Package Size: #AP80311-1 50ug #AP80311-2 100ug #AP80311-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	ADAM 10,Kuzbanian protein homolog,Mammalian disintegrin-metalloprotease,CD156c
Accession No.	Uniprot ID:O77633
Uniprot	O77633
GeneID	397345;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADAM10

Note: This product is for in vitro research use only