

Activity-dependent neuroprotector homeobox protein

Catalog No: #AP80313



Package Size: #AP80313-1 50ug #AP80313-2 100ug #AP80313-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Activity-dependent neuroprotector homeobox protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Activity-dependent neuroprotective protein,ADNP1,KIAA0784
Accession No.	Uniprot ID:Q9H2P0
Uniprot	Q9H2P0
GeneID	23394;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADNP

Note: This product is for in vitro research use only