

# Disintegrin and metalloproteinase domain-containing protein 17

Catalog No: #AP80315

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP80315-1 50ug #AP80315-2 100ug #AP80315-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ADAM 17, TNF-alpha convertase, TNF-alpha-converting enzyme, CD156b, Tace
Accession No.	Uniprot ID: Q9Z1K9
Uniprot	Q9Z1K9
GeneID	57027;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Adam17

Note: This product is for in vitro research use only