

# Disintegrin and metalloproteinase domain-containing protein 9



Catalog No: #AP80316

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP80316-1 50ug #AP80316-2 100ug #AP80316-3 1mg

## Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADAM 9, Meltrin-gamma, Metalloprotease/disintegrin/cysteine-rich protein 9, Myeloma cell metalloproteinase, Kiaa0021, Mdc9, Mltng
Accession No.	Uniprot ID:Q61072
Uniprot	Q61072
GeneID	11502;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Adam9

Note: This product is for in vitro research use only