

Beta-1 adrenergic receptor

Catalog No: #AP80318



Package Size: #AP80318-1 50ug #AP80318-2 100ug #AP80318-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-1 adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Beta-1 adrenoreceptor,Beta-1 adrenoceptor,Adrb1r
Accession No.	Uniprot ID:P34971
Uniprot	P34971
GeneID	11554;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adrb1

Note: This product is for in vitro research use only