

Beta-2 adrenergic receptor

Catalog No: #AP80322



Package Size: #AP80322-1 50ug #AP80322-2 100ug #AP80322-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beta-2 adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Beta-2 adrenoreceptor,Beta-2 adrenoceptor,ADRB2R,B2AR
Accession No.	Uniprot ID:P07550
Uniprot	P07550
GeneID	154;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADRB2

Note: This product is for in vitro research use only