

Delta-1-pyrroline-5-carboxylate dehydrogenase, mitochondrial

Catalog No: #AP80323

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80323-1 50ug #AP80323-2 100ug #AP80323-3 1mg

Description

| | |
|--------------------|--|
| Product Name | Delta-1-pyrroline-5-carboxylate dehydrogenase, mitochondrial |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | P5C dehydrogenase,Aldehyde dehydrogenase family 4 member A1,L-glutamate gamma-semialdehyde dehydrogenase |
| Accession No. | Uniprot ID:Q8CHT0 |
| Uniprot | Q8CHT0 |
| GeneID | 212647; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aldh4a1

Note: This product is for in vitro research use only