

Aminoacylase-1

Catalog No: #AP80324



Package Size: #AP80324-1 50ug #AP80324-2 100ug #AP80324-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aminoacylase-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ACY-1,N-acyl-L-amino-acid amidohydrolase
Accession No.	Uniprot ID:Q03154
Uniprot	Q03154
GeneID	95;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ACY1

Note: This product is for in vitro research use only