

Adenosine deaminase

Catalog No: #AP80326



Package Size: #AP80326-1 50ug #AP80326-2 100ug #AP80326-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adenosine deaminase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Adenosine aminohydrolase
Accession No.	Uniprot ID:Q920P6
Uniprot	Q920P6
GeneID	24165;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ada

Note: This product is for in vitro research use only