

Retinal dehydrogenase 1

Catalog No: #AP80327



Package Size: #AP80327-1 50ug #AP80327-2 100ug #AP80327-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Retinal dehydrogenase 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | RALDH 1,ALDH-E1,ALHDII,Aldehyde dehydrogenase family 1 member A1,Aldehyde dehydrogenase, cytosolic,Ahd-2,Ahd2,Aldh1 |
| Accession No. | Uniprot ID:P24549 |
| Uniprot | P24549 |
| GeneID | 11668; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aldh1a1

Note: This product is for in vitro research use only