

Serum albumin

Catalog No: #AP80328



Package Size: #AP80328-1 50ug #AP80328-2 100ug #AP80328-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serum albumin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GIG20,GIG42,PRO0903,PRO1708,PRO2044,PRO2619,PRO2675,UNQ696/PRO1341
Accession No.	Uniprot ID:P02768
Uniprot	P02768
GeneID	213;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALB

Note: This product is for in vitro research use only