

ATP-citrate synthase

Catalog No: #AP80329

Package Size: #AP80329-1 50ug #AP80329-2 100ug #AP80329-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP-citrate synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ATP-citrate (pro-S-)-lyase,ACL,Citrate cleavage enzyme
Accession No.	Uniprot ID:P53396
Uniprot	P53396
GeneID	47;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ACLY

Note: This product is for in vitro research use only