

Long-chain-fatty-acid--CoA ligase 1

Catalog No: #AP80333

Package Size: #AP80333-1 50ug #AP80333-2 100ug #AP80333-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Long-chain-fatty-acid--CoA ligase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Long-chain acyl-CoA synthetase 1,LACS 1,Long-chain-fatty-acid--CoA ligase, liver isozyme,Acs1,Acsl2,Facl2
Accession No.	Uniprot ID:P18163
Uniprot	P18163
GeneID	25288;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Acsl1

Note: This product is for in vitro research use only