

Alpha-2A adrenergic receptor

Catalog No: #AP80334

Package Size: #AP80334-1 50ug #AP80334-2 100ug #AP80334-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-2A adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-2A adrenoreceptor,Alpha-2A adrenoceptor
Accession No.	Uniprot ID:Q01338
Uniprot	Q01338
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adra2a

Note: This product is for in vitro research use only