

Serine/threonine-protein kinase receptor R3

Catalog No: #AP80336



Package Size: #AP80336-1 50ug #AP80336-2 100ug #AP80336-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase receptor R3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SKR3,Activin receptor-like kinase 1,ALK-1,TGF-B superfamily receptor type I,TSR-I,Acvrlk1,Alk-1
Accession No.	Uniprot ID:Q61288
Uniprot	Q61288
GeneID	11482;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Acvrl1

Note: This product is for in vitro research use only