

# Disintegrin and metalloproteinase domain-containing protein 29



Catalog No: #AP80343

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP80343-1 50ug #AP80343-2 100ug #AP80343-3 1mg

## Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 29
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADAM 29
Accession No.	Uniprot ID:Q811Q4
Uniprot	Q811Q4
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Adam29

Note: This product is for in vitro research use only