

## Protein argonaute-2

Catalog No: #AP80344



Package Size: #AP80344-1 50ug #AP80344-2 100ug #AP80344-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein argonaute-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Argonaute2,Argonaute RISC catalytic component 2,Eukaryotic translation initiation factor 2C 2,eIF-2C 2,PAZ Piwi domain protein,PPD,Protein slicer,EIF2C2
Accession No.	Uniprot ID:Q9UKV8
Uniprot	Q9UKV8
GeneID	27161;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:AGO2

Note: This product is for in vitro research use only