

# Type-1 angiotensin II receptor

Catalog No: #AP80345



Package Size: #AP80345-1 50ug #AP80345-2 100ug #AP80345-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Type-1 angiotensin II receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Angiotensin II type-1 receptor,AT1
Accession No.	Uniprot ID:P30555
Uniprot	P30555
GeneID	100519976;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:AGTR1

Note: This product is for in vitro research use only