

Agrin

Catalog No: #AP80347



Package Size: #AP80347-1 50ug #AP80347-2 100ug #AP80347-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Agrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Agrin
Accession No.	Uniprot ID:A2ASQ1
Uniprot	A2ASQ1
GeneID	11603;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Agrin

Note: This product is for in vitro research use only