

# Rubisco (Large Chain) Monoclonal Antibody(1H1)

Catalog No: #44092

Package Size: #44092-1 50ul #44092-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

|                       |   |
|-----------------------|---|
| Product Name          | Rubisco (Large Chain) Monoclonal Antibody(1H1)  |
| Host Species          | Mouse   |
| Clonality             | Monoclonal  |
| Clone No.             | 1H1   |
| Isotype               | IgG1  |
| Applications          | WB  |
| Species Reactivity    | Plant   |
| Specificity           | Rubisco(Large Chain) Mouse Monoclonal Antibody(1H1) detects endogenous levels of Rubisco(Large Chain) |
| Immunogen Type        | Recombinant protein   |
| Immunogen Description | Recombinant Protein   |
| Target Name           | Rubisco (Large Chain)   |
| Calculated MW         | 53kd  |
| Concentration         | 1mg/ml  |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                               |
| Storage               | Store at -20°C/1 year   |

## Application Details

Western blotting: 1:2000-1:5000

## Images



Western blot analysis of Arabidopsis with Rubisco(Large Chain) Mouse mAb diluted at  
 1)1:2,000  
 2)1:5,000

Note: This product is for in vitro research use only