

Adenylate cyclase type 8

Catalog No: #AP80353



Package Size: #AP80353-1 50ug #AP80353-2 100ug #AP80353-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenylate cyclase type 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ATP pyrophosphate-lyase 8, Adenylate cyclase type VIII, Adenylyl cyclase 8, Ca(2+)/calmodulin-activated adenylyl cyclase
Accession No.	Uniprot ID:P97490
Uniprot	P97490
GeneID	11514;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Adcy8

Note: This product is for in vitro research use only