

Aldo-keto reductase family 1 member C4

Catalog No: #AP80355



Package Size: #AP80355-1 50ug #AP80355-2 100ug #AP80355-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aldo-keto reductase family 1 member C4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	3-alpha-HSD1,3-alpha-hydroxysteroid dehydrogenase type I,Chlordecone reductase,CDR,Dihydrodiol dehydrogenase 4,DD-4,HAKRA,CHDR
Accession No.	Uniprot ID:P17516
Uniprot	P17516
GeneID	1109;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AKR1C4

Note: This product is for in vitro research use only