

RAC-beta serine/threonine-protein kinase

Catalog No: #AP80356



Package Size: #AP80356-1 50ug #AP80356-2 100ug #AP80356-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	RAC-beta serine/threonine-protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein kinase Akt-2,Protein kinase B beta,PKB beta,RAC-PK-beta
Accession No.	Uniprot ID:Q60823
Uniprot	Q60823
GeneID	11652;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Akt2

Note: This product is for in vitro research use only