

Fructose-bisphosphate aldolase C

Catalog No: #AP80358



Package Size: #AP80358-1 50ug #AP80358-2 100ug #AP80358-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fructose-bisphosphate aldolase C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Brain-type aldolase,ALDC
Accession No.	Uniprot ID:P09972
Uniprot	P09972
GeneID	230;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALDOC

Note: This product is for in vitro research use only