

## TBP/TATA Binding Protein Monoclonal Antibody(4H2)

Catalog No: #44093

Package Size: #44093-1 50ul #44093-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TBP/TATA Binding Protein Monoclonal Antibody(4H2)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4H2
Isotype	IgG1
Applications	WB
Species Reactivity	Hu Rt Ms
Specificity	TBP/TATA Binding Protein Mouse Monoclonal Antibody(4H2) detects endogenous levels of TBP/TATA Binding Protein
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	TBP/TATA Binding Protein
Accession No.	Swiss-Prot:P20226Gene ID:6908
Uniprot	P20226
GeneID	6908;
Calculated MW	35-45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:1000-1:2000

## Images



Western blot analysis of 1) Hela, 2) Mouse Brain Tissue, 3) Rat Brain Tissue with TBP/TATA Binding Protein Mouse mAb diluted at 1:2,000.

---

Note: This product is for in vitro research use only