

# Disintegrin and metalloproteinase domain-containing protein 10



Catalog No: #AP80363

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP80363-1 50ug #AP80363-2 100ug #AP80363-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADAM 10,Kuzbanian protein homolog,Mammalian disintegrin-metalloprotease,CD156c,Kuz,Madm
Accession No.	Uniprot ID:O35598
Uniprot	O35598
GeneID	11487;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Adam10

Note: This product is for in vitro research use only