

ADP/ATP translocase 3

Catalog No: #AP80365



Package Size: #AP80365-1 50ug #AP80365-2 100ug #AP80365-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADP/ATP translocase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADP,ATP carrier protein 3,ADP,ATP carrier protein, isoform T2,ANT 2,Adenine nucleotide translocator 3,ANT 3,Solute carrier family 25 member 6,ANT3,CDABP0051
Accession No.	Uniprot ID:P12236
Uniprot	P12236
GeneID	293;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC25A6

Note: This product is for in vitro research use only