

## Acyl-CoA desaturase

Catalog No: #AP80368



Package Size: #AP80368-1 50ug #AP80368-2 100ug #AP80368-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Acyl-CoA desaturase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Delta(9)-desaturase,Delta-9 desaturase,Fatty acid desaturase,Stearoyl-CoA desaturase,hSCD1
Accession No.	Uniprot ID:O00767
Uniprot	O00767
GeneID	6319;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SCD

Note: This product is for in vitro research use only