

Acrosin

Catalog No: #AP80369



Package Size: #AP80369-1 50ug #AP80369-2 100ug #AP80369-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Acrosin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	53 kDa fucose-binding protein
Accession No.	Uniprot ID:P08001
Uniprot	P08001
GeneID	397098;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ACR

Note: This product is for in vitro research use only