

## Cyclophilin B Monoclonal Antibody(2B10)

Catalog No: #44094

Package Size: #44094-1 50ul #44094-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Cyclophilin B Monoclonal Antibody(2B10)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2B10
Isotype	IgG1
Applications	WB
Species Reactivity	Hu Rt Ms
Specificity	Cyclophilin B Mouse Monoclonal Antibody(2B10) detects endogenous levels of Cyclophilin B
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Cyclophilin B
Accession No.	Swiss-Prot:P23284Gene ID:5479
Uniprot	P23284
GeneID	5479;
Calculated MW	21kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:1000-1:2000

## Images



Western blot analysis of 1) Jurkat, 2)293T, 3)Rat Liver Tissue, 4)3T3 with Cyclophilin B Mouse mAb diluted at 1:2,000.

---

Note: This product is for in vitro research use only