

# 5'-AMP-activated protein kinase catalytic subunit alpha-2

Catalog No: #AP80372

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP80372-1 50ug #AP80372-2 100ug #AP80372-3 1mg

## Description

Product Name	5'-AMP-activated protein kinase catalytic subunit alpha-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AMPK subunit alpha-2,Acetyl-CoA carboxylase kinase,ACACA kinase,Hydroxymethylglutaryl-CoA reductase kinase,HMGCR kinase
Accession No.	Uniprot ID:Q8BRK8
Uniprot	Q8BRK8
GeneID	108079;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Prkaa2

Note: This product is for in vitro research use only