

ATP-binding cassette sub-family B member 9

Catalog No: #AP80373



Package Size: #AP80373-1 50ug #AP80373-2 100ug #AP80373-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP-binding cassette sub-family B member 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ATP-binding cassette transporter 9,ABC transporter 9 protein,TAP-like protein,TAPL
Accession No.	Uniprot ID:Q9QYJ4
Uniprot	Q9QYJ4
GeneID	63886;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abcb9

Note: This product is for in vitro research use only