

# Beta-2-glycoprotein 1

Catalog No: #AP80379

Package Size: #AP80379-1 50ug #AP80379-2 100ug #AP80379-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| Product Name       | Beta-2-glycoprotein 1                                  |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli                                                 |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Bovine                                                 |
| Other Names        | Apolipoprotein H,Apo-H,Beta-2-glycoprotein I,B2GPI     |
| Accession No.      | Uniprot ID:P17690                                      |
| Uniprot            | P17690                                                 |
| GeneID             | 281006;                                                |
| Tag Info           | His                                                    |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:APOH

Note: This product is for in vitro research use only