

Desmin Monoclonal Antibody(1B12)

Catalog No: #44095

Package Size: #44095-1 50ul #44095-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | Desmin Monoclonal Antibody(1B12) |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 1B12 |
| Isotype | IgG1 |
| Applications | IHC |
| Species Reactivity | Hu Rt Ms |
| Specificity | Desmin Mouse Monoclonal Antibody(1B12) detects endogenous levels of Desmin |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Recombinant Protein |
| Target Name | Desmin |
| Accession No. | Swiss-Prot:P17661Gene ID:1674 |
| Uniprot | P17661 |
| GeneID | 1674; |
| Calculated MW | 53kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Mouse Skeletal Muscle Tissue using DesminB&B™mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Heart Tissue using DesminMouse mAb diluted at 1:200.

Note: This product is for in vitro research use only