

Alcohol dehydrogenase 1

Catalog No: #AP80380



Package Size: #AP80380-1 50ug #AP80380-2 100ug #AP80380-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alcohol dehydrogenase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADH-A2,Alcohol dehydrogenase A subunit,Adh-1
Accession No.	Uniprot ID:P00329
Uniprot	P00329
GeneID	11522;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adh1

Note: This product is for in vitro research use only