

Activity-dependent neuroprotector homeobox protein

Catalog No: #AP80381



Package Size: #AP80381-1 50ug #AP80381-2 100ug #AP80381-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Activity-dependent neuroprotector homeobox protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Activity-dependent neuroprotective protein
Accession No.	Uniprot ID:Q9JKL8
Uniprot	Q9JKL8
GeneID	64622;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Adnp

Note: This product is for in vitro research use only