

Protein argonaute-2

Catalog No: #AP80384



Package Size: #AP80384-1 50ug #AP80384-2 100ug #AP80384-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein argonaute-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Argonaute2,Argonaute RISC catalytic component 2,Eukaryotic translation initiation factor 2C 2,eIF-2C 2,Piwi/argonaute family protein melf2c2,Protein slicer,Eif2c2,Kiaa4215
Accession No.	Uniprot ID:Q8CJG0
Uniprot	Q8CJG0
GeneID	239528;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ago2

Note: This product is for in vitro research use only