

Activation-induced cytidine deaminase

Catalog No: #AP80386



Package Size: #AP80386-1 50ug #AP80386-2 100ug #AP80386-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Activation-induced cytidine deaminase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Cytidine aminohydrolase,AID,Cytidine aminohydrolase
Accession No.	Uniprot ID:Q75W64
Uniprot	Q75W64
GeneID	442983;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AICDA

Note: This product is for in vitro research use only