

Angiopoietin-4

Catalog No: #AP80388



Package Size: #AP80388-1 50ug #AP80388-2 100ug #AP80388-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Angiopoietin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ANG-4,Angiopoietin-3,ANG-3,Agpt4,Ang3
Accession No.	Uniprot ID:Q9WVH6
Uniprot	Q9WVH6
GeneID	11602;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Angpt4

Note: This product is for in vitro research use only