

Amelogenin, X isoform

Catalog No: #AP80393

Package Size: #AP80393-1 50ug #AP80393-2 100ug #AP80393-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Amelogenin, X isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Class I amelogenin
Accession No.	Uniprot ID:P02817
Uniprot	P02817
GeneID	281620;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: AMELX

Note: This product is for in vitro research use only