

AKT-interacting protein

Catalog No: #AP80395

Package Size: #AP80395-1 50ug #AP80395-2 100ug #AP80395-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AKT-interacting protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ft1,Fused toes protein homolog,FTS
Accession No.	Uniprot ID:Q9H8T0
Uniprot	Q9H8T0
GeneID	64400;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AKTIP

Note: This product is for in vitro research use only