

Delta-1-pyrroline-5-carboxylate dehydrogenase, mitochondrial

Catalog No: #AP80396

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80396-1 50ug #AP80396-2 100ug #AP80396-3 1mg

Description

Product Name	Delta-1-pyrroline-5-carboxylate dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	P5C dehydrogenase,Aldehyde dehydrogenase family 4 member A1,L-glutamate gamma-semialdehyde dehydrogenase
Accession No.	Uniprot ID:P0C2X9
Uniprot	P0C2X9
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aldh4a1

Note: This product is for in vitro research use only