

Disintegrin and metalloproteinase domain-containing protein 10



Catalog No: #AP80397

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80397-1 50ug #AP80397-2 100ug #AP80397-3 1mg

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ADAM 10,Kuzbanian protein homolog,Mammalian disintegrin-metalloprotease,Myelin-associated metalloproteinase,CD156c,MADM
Accession No.	Uniprot ID:Q10741
Uniprot	Q10741
GeneID	282132;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADAM10

Note: This product is for in vitro research use only