

Adiponectin

Catalog No: #AP80398



Package Size: #AP80398-1 50ug #AP80398-2 100ug #AP80398-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adiponectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Adipocyte-specific protein AdipoQ, Adipocyte, C1q and collagen domain-containing protein, 30 kDa adipocyte complement-related protein, Adipocyte complement-related 30 kDa protein (ACRP30)
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ADIPOQ

Note: This product is for in vitro research use only