

Apolipoprotein L1

Catalog No: #AP80400



Package Size: #AP80400-1 50ug #AP80400-2 100ug #AP80400-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein L1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Apolipoprotein L, Apo-L, Apolipoprotein L-I, ApoL-I, APOL
Accession No.	Uniprot ID:O14791
Uniprot	O14791
GeneID	8542;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: APOL1

Note: This product is for in vitro research use only